

FINALLY! Non-Cycling Alarm for the VAR™

VAR™ - Monitor

- Works with current VAR models
- Audible and visible alarm
- Operates with a 9 VDC battery
- Easy to setup and use
- Cost effective for emergency preparedness



Order No. VM-3500



The **VAR™-Monitor** is a battery operated device that can be attached to any VAR™ (VORTAN Automatic Resuscitator) to monitor cycling conditions during operation. The device will check for the non-cycling condition of the VAR™ and activate both visual and audible alarms if the pressure remains unchanged for more than the preset time of 8 seconds. The **VAR™-Monitor** is ideal for managing multiple patients on VAR™ in any ventilation emergency.



21 Goldenland Court, #100 • Sacramento, CA 95834
(800) 434-4034 • tx (916) 648-9751 • www.vortan.com