

Repair not Replace

when a pilot line becomes damaged!

Our **Pilot Tube Repair Kit**, allows you to repair a damaged pilot line on an endotracheal or tracheotomy tube.

Fits pilot line of most tubes from **0.032 in. O.D. (21g)** to **0.050 in. O.D. (18g)**.

Can also be used with inflatable face masks.

Essential for all critical care areas.

*Clamp,
cut, insert...*

*Unclamp,
reinflate...*

That's it!



**Instrumentation
Industries, Inc.**

1-800-633-8577

Order Online! www.iiimedical.com



ISO 13485: 2003
FM 72200